



IF IT WORKS, IT IS NOT STUPID.

BUSTA LAB MOTTO 1

JUST BECAUSE SOMETHING HAS ALWAYS BEEN DONE ONE WAY DOES NOT MEAN IT  
CAN'T BE IMPROVED.

BUSTA LAB MOTTO 2

THE BUSTA LAB

# THE BUSTA LAB HANDBOOK

Copyright © 2021 The Busta Lab

PUBLISHED BY

TUFTE-LATEX.GITHUB.IO/TUFTE-LATEX/

Licensed under the Apache License, Version 2.0 (the “License”); you may not use this file except in compliance with the License. You may obtain a copy of the License at <http://www.apache.org/licenses/LICENSE-2.0>. Unless required by applicable law or agreed to in writing, software distributed under the License is distributed on an “AS IS” BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied. See the License for the specific language governing permissions and limitations under the License.

*First printing, December 2021*

# *Contents*

*Nucleic Acids*      13

*Materials and Reagents*      17

*Bibliography*      19



## *List of Figures*





## *List of Tables*



*Dedicated to those who appreciate  $\text{\LaTeX}$  and the work of  
Edward R. Tufte and Donald E. Knuth.*



## *Nucleic Acids*

## *Tissue Harvesting and Homogenization*

### **Tissue harvesting and homogenization**

- [ ] Isolate tissue with RNase-free tweezers and scissors, put in 50 ml falcon tube, flash freeze.
- [ ] Transfer small aliquots of tissue from 50 ml falcon tube to mortar, grind to a *FINE* powder, liquid N<sub>2</sub> as needed.
- [ ] Use pre-labelled pre-frozen 2 ml tubes to scoop and store powder, always keeping tissue frozen. Store at -80°C.

[illegible]

## mRNA Sequencing

### cDNA library preparation (7 hours)

#### Reverse Transcription

[ ] *Mix 1*: 1 ng mRNA (x  $\mu$ l), 1  $\mu$ l **VNP**, 1  $\mu$ l dNTPs, 9-x  $\mu$ l RNase-free water, FMSD<sup>6</sup>, 65°C 5 min, snap freeze.<sup>7</sup>

[ ] *Mix 2*: 4  $\mu$ l RT buffer, 1  $\mu$ l RNaseOUT, 1  $\mu$ l RNase-free water, 2  $\mu$ l **SSP**, FMSD.

[ ] Make *RT Mix 1*: Add *Mix 2* to *Mix 1*, FMSD, Incubate (42°C, 3 min).

[ ] Make *RT Mix 2*: Add 1  $\mu$ l RT enzyme<sup>8</sup> (200 U/ $\mu$ l) to *RT Mix 1*. FMSD. Incubate (42°C, 90 min).

[ ] Inactivate *RT Mix 2* by incubating (85°C, 5 min; hold at 4°C).

\*\* Is there some way to run QC on length distribution here???

I abandoned cDNA kit because I was getting an abundance of short reads 150bp to 300bp that I didn't even see with direct RNA seq which doesn't have a strand switch oligo to select full length transcripts'.

I ended up switching to the direct RNA kit, which confirmed it wasn't the input RNA's fault as the read length counts distribution looked great compared to direct cDNA. And the distribution was correlated well with what I saw on the bioanalyzer. Just in general Nanopore need to do more on the RNA seq side of things.

Reducing PCR cycles (if possible) can be useful, PCR tends to favour shorter amplicons as they amplify more efficiently. So doing as few cycles as possible can help give you a higher proportion of longer reads.

#### PCR amplification (2 hours)

[ ] Make 4 identical tubes of *Mix 3*: 18.5  $\mu$ l nuc-free water, 25  $\mu$ l 2x LongAmp MM, 1.5  $\mu$ l **cPRM**, 5  $\mu$ l *RT Mix 2*. FMSD.

[ ] Run PCR on *Mix 3*: [95°C - 1 min], [95°C, 15s | 62°C, 15s | 65°C, 150s] x 30 cycles<sup>9</sup>, [65°C, 6 min], [4°C hold]. \$

[ ] Add 1  $\mu$ l Exonuclease (20 U) to each tube of *Mix 3*, FMSD, 37°C - 15 min, 80°C - 15 min (degrades primers).

#### Cleanup and Concentrate cDNA (1 hour)

[ ] Resuspend beads by vortexing for 1 min and let them warm to room temperature.

[ ] Pool *Mix 3s* in new 1.5ml LoBind. Vortex beads (10 s), add 120  $\mu$ l beads to *Mix 3*.<sup>10</sup> Rock (5 min, RT).

[ ] Prepare 500  $\mu$ l fresh 70% ethanol using RNase-free water (350  $\mu$ l EtOH + 150  $\mu$ l water).

[ ] Pellet beads on magnet at least 5 minutes [DO NOT spin down], remove supernatant. Prepare gel for quality check.

[ ] [ ] [ ] 3x: Tube on magnet, add 1000 (/enough to cover pellet)  $\mu$ l 70% EtOH, beads sit in EtOH for 30s, discard EtOH in two batches.

[ ] Tube on magnet: with a P10, remove residual EtOH, then dry for 30 seconds.

[ ] Remove tube from rack, resuspend gently by pipetting in 12  $\mu$ l Elution Buffer. Incubate RT 10 minutes.

[ ] Pellet beads on magnet 5 min [DO NOT spin down] remove + retain 12  $\mu$ l of eluate in new 1.5 mL LoBind, spent beads to biohazard.

[ ] Analyze 1  $\mu$ l of amplified cDNA for size, quantity, quality using Qubit and gel: ng/ $\mu$ l,  $\mu$ l total, and bp avg size. \$

### Sequencing

#### Rapid Adapter Addition

[ ] Thaw Sequencing buffer (SQB), Loading Beads (LB), Flush tether (FLT), and one tube flush buffer (PFB/FB).

[ ] Make up 100-200 fmol cDNA in Elution Buffer (EB; final vol. 12  $\mu$ l). At 1kb, 100 fmol = 60 ng, 200 fmol = 120 ng.

[ ] Add Rapid Adapter (RAP; 1  $\mu$ l), mix by pipetting, spin down. Incubate (10+ min RT). 4.2 On MinION Flowcell Run on Nanopore flowcell (1 hr + up to 48 hours runtime)

#### Sequencing on MinION

[ ] Set P1000 to 200  $\mu$ l and insert tip into priming port. Using pipette suction method, turn the dial to 220-230  $\mu$ l, or until you can see a small amount of buffer entering the pipette tip. [ ] Prepare the flow cell priming mix: 30  $\mu$ l of thawed and mixed Flush tether (FLT) directly to the thawed and mixed Flush buffer (FB). Mix by vortexing. [ ] With SpotON closed, add 800  $\mu$ l PriMix to the priming port without introducing air into the system. Wait 10 minutes. [ ]

Add 37.5  $\mu$ l SQB to the cDNA, thoroughly mix Loading Beads by pipetting, add 25.5  $\mu$ l beads to 12  $\mu$ l cDNA to make Sample. [ ] Add 200  $\mu$ l PriMix to the priming port without introducing air into the system. [ ] Mix the prepared library gently by pipetting just prior to loading. [ ] Open SpotON, add 75  $\mu$ l Sample dropwise to the SpotON port. Ensure each drop flows into the cell before adding the next. Close SpotON, close priming port, close lid.

*Sequencing on Flongle* To load the Flongle flow cell, 7.5  $\mu$ L of sequencing buffer (Oxford Nanopore Technologies), 5  $\mu$ L of freshly mixed loading beads (Oxford Nanopore Technologies), and 2.5  $\mu$ L of the prepped DNA were mixed in a new 1.5 mL Eppendorf tube.

Following priming, the freshly mixed solution of sequencing buffer, loading beads, and DNA library was added in a dropwise fashion to the sample port of the flow cell, and a new sequencing run was started. If pore occupancy was initially low (i.e. < 10 pores sequencing), a second library (again consisting of 7.5  $\mu$ L sequencing buffer, 5  $\mu$ L loading beads, and 2.5  $\mu$ L DNA) was added to increase sequencing yield.

<sup>6</sup> FMSD = flick mix and spin down

<sup>7</sup> on block

<sup>8</sup> Superscript Reverse Transcriptase IV or Maxima H Minus Reverse Transcriptase

<sup>9</sup> Nanopore says 11-22 cycles, start with 14. We find 30 is usually necessary, though too many cycles can create "PCR duplicates" apparently

<sup>10</sup> Bead:DNA ratio is critical!! 120  $\mu$ l of beads is a 0.6 ratio which should give 400-500bp+ in the eluate. 100  $\mu$ l of beads is a 0.5 ratio which should give 600bp+ in the eluate. Other ratios are also possible.



# Materials and Reagents

## Reagents

**dNTPs**  
(10mM)

## Nanopore Primers

### VNP

(2uM) 5' - 5phosACTTGCCTGTCGCTCTATCTTCTTTTTTTTTTTTTTTTTTVN - 3' Where V = A, C, or G, and N = A, C, G, or T The 5' phosphate is added as a modification during oligo synthesis for improved ligation efficiency. However, this phosphate is not essential for the Direct cDNA Sequencing protocol, since the end-repair step also adds a 5' phosphate.

### SSP

5' - TTTCTGTTGGTGCTGATATTGCT mGmGmG - 3'

### cPRM

Forward: 5' - atcgctaccgtgacaagaaagttgctggtgtctttgtgACTTGCCTGTCGCTCTATCTTC - 3' 5' - atcgctaccgtgacaagaaagttgctggtgtctttgtgAGCTAGcACTTGCCTGTCGCTCTATCTTC - 3' Reverse: 5' - atcgctaccgtgacaagaaagttgctggtgtctttgtgTTTCTGTTGGTGCTGATATTGC - 3' 5' - atcgctaccgtgacaagaaagttgctggtgtctttgtgAGCTAGcTTTCTGTTGGTGCTGATATTGC - 3'

TO MAKE CUSTOM BARCODES, INSERT CUSTOM SEQUENCE BETWEEN CAPS/LOWER CASE ABOVE



## *Bibliography*

- [1] Piotr Chomczynski and Nicoletta Sacchi, "The single-step method of RNA isolation by acid guanidinium thiocyanate-phenol-chloroform extraction: twenty-something years on", *Nature Protocols*, 1(2), 2006, pp. 581–585.